



Sample Name: CZ126S-30_HexIPA_010_AD

```
=====
Acq. Operator   : AK Bräse                      Seq. Line :   11
Acq. Instrument : Instrument 1                  Location  : Vial 14
Injection Date  : 12/22/2019 6:25:34 AM        Inj       :    1
                                                Inj Volume: 15 µl
Different Inj Volume from Sequence !      Actual Inj Volume : 1 µl
Acq. Method     : C:\CHEM32\1\METHODS\METHODS.OLD\CZ_010_20MIN_ISO_CHIRAL.M
Last changed    : 12/7/2019 1:10:40 AM by AK Bräse
Analysis Method : C:\CHEM32\1\METHODS\YMC_CZ_HEX-IPR_100-0_20MIN_1ML_ISOCHIR.M
Last changed    : 12/8/2019 10:11:59 PM by AK Bräse
                  (modified after loading)
Method Info     : YMC Amylose-SA, 4.6 x 250 mm, 5 µm,
                  n-hexane, 1 mL/min
=====
```



Sample Name: CZ126S-30_HexIPA_010_AD

=====
Area Percent Report
=====

Sorted By : Signal
Multiplier: : 1.0000
Dilution: : 1.0000
Use Multiplier & Dilution Factor with ISTDs

Sample Name: CZ126S-30_HexIPA_010_AD

Signal 1: DAD1 A, Sig=200,4 Ref=400,4

Peak #	RetTime [min]	Type	Width [min]	Area [mAU*s]	Height [mAU]	Area %
1	4.216	BV	0.0869	1774.04712	318.79095	50.8888
2	4.513	VB	0.0895	1712.07471	295.74234	49.1112

Totals : 3486.12183 614.53329

Signal 2: DAD1 B, Sig=256,4 Ref=400,4

Signal 3: DAD1 C, Sig=330,4 Ref=400,4

Signal 4: DAD1 D, Sig=450,4 Ref=600,4

Signal 5: DAD1 E, Sig=560,4 Ref=600,4

=====
*** End of Report ***